

# Recombinant Human TLE1

Catalog No: #GP13052



Package Size: #GP13052-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TLE1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 571-770 amino acids of human TLE1
Target Name	transducin like enhancer of split 1
Other Names	ESG; ESG1; GRG1
Accession No.	Swissprot:Q04724Gene Accession:BC015747
Uniprot	Q04724
GeneID	7088;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only