

Recombinant Human MIEN1

Catalog No: #GP13070



Package Size: #GP13070-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MIEN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-112 amino acids of human MIEN1
Target Name	migration and invasion enhancer 1
Other Names	C35; ORB3; XTP4; RDX12; C17orf37
Accession No.	Swissprot:Q9BRT3Gene Accession:BC006006
Uniprot	Q9BRT3
GeneID	84299;
Storage	-20~-80°C, pH 7.6 PBS

Background

Increases cell migration by inducing filopodia formation at the leading edge of migrating cells. Plays a role in regulation of apoptosis, possibly through control of CASP3. May be involved in a redox-related process.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only