

# Recombinant Human AP1S1

Catalog No: #GP13078



Package Size: #GP13078-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human AP1S1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-133 amino acids of human AP1S1
Target Name	adaptor related protein complex 1 sigma 1 subunit
Other Names	AP19; EKV3; CLAPS1; MEDNIK; SIGMA1A
Accession No.	Swissprot:P61966Gene Accession:BC003561
Uniprot	P61966
GeneID	1174;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is part of the clathrin coat assembly complex which links clathrin to receptors in coated vesicles. These vesicles are involved in endocytosis and Golgi processing. This protein, as well as beta-prime-adaptin, gamma-adaptin, and the medium (mu) chain AP47, form the AP-1 assembly protein complex located at the Golgi vesicle.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only