

Recombinant Human PTP4A3

Catalog No: #GP13083



Package Size: #GP13083-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PTP4A3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-148 amino acids of human PTP4A3
Target Name	protein tyrosine phosphatase type IVA, member 3
Other Names	PRL3; PRL-3; PRL-R
Accession No.	Swissprot:O75365Gene Accession:BC003105
Uniprot	O75365
GeneID	11156;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the protein-tyrosine phosphatase family. Protein tyrosine phosphatases are cell signaling molecules that play regulatory roles in a variety of cellular processes. Studies of this class of protein tyrosine phosphatase in mice demonstrates that they are prenylated in vivo, suggesting their association with cell plasma membrane. The encoded protein may enhance cell proliferation, and overexpression of this gene has been implicated in tumor metastasis. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only