

Recombinant Human MRPL20

Catalog No: #GP13085



Package Size: #GP13085-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MRPL20
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 46-149 amino acids of human MRPL20
Target Name	mitochondrial ribosomal protein L20
Other Names	L20mt; MRP-L20
Accession No.	Swissprot:Q9BYC9Gene Accession:BC009515
Uniprot	Q9BYC9
GeneID	55052;
Storage	-20~-80°C, pH 7.6 PBS

Background

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. A pseudogene corresponding to this gene is found on chromosome 21q. Alternative splicing results in multiple transcript variants encoding different isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only