

Recombinant Human THAP12

Catalog No: #GP13086



Package Size: #GP13086-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human THAP12
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 612-761 amino acids of human THAP12
Target Name	THAP domain containing 12
Other Names	DAP4; THAP0; PRKRIR; P52rIPK
Accession No.	Swissprot:O43422Gene Accession:BC021992
Uniprot	O43422
GeneID	5612;
Storage	-20~-80°C, pH 7.6 PBS

Background

Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only