

# Recombinant Human C8orf44

Catalog No: #GP13092



Package Size: #GP13092-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human C8orf44
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-159 amino acids of human C8orf44
Target Name	chromosome 8 open reading frame 44
Accession No.	Swissprot:Q96CB5Gene Accession:BC014448
Uniprot	Q96CB5
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only