

Recombinant Human DGLUCY

Catalog No: #GP13093



Package Size: #GP13093-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DGLUCY
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-159 amino acids of human C14orf159
Target Name	chromosome 14 open reading frame 159
Accession No.	Swissprot:Q7Z3D6Gene Accession:BC009182
Uniprot	Q7Z3D6
GeneID	80017;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only