

Recombinant Human ISG15

Catalog No: #GP13095



Package Size: #GP13095-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ISG15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-157 amino acids of human ISG15
Target Name	ISG15 ubiquitin-like modifier
Other Names	G1P2; IP17; UCRP; IFI15; IMD38; hUCRP
Accession No.	Swissprot:P05161Gene Accession:BC009507
Uniprot	P05161
GeneID	9636;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a ubiquitin-like protein that is conjugated to intracellular target proteins upon activation by interferon-alpha and interferon-beta. Several functions have been ascribed to the encoded protein, including chemotactic activity towards neutrophils, direction of ligated target proteins to intermediate filaments, cell-to-cell signaling, and antiviral activity during viral infections. While conjugates of this protein have been found to be noncovalently attached to intermediate filaments, this protein is sometimes secreted.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only