

TBC1D1 Antibody

Catalog No: #43157

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBC1D1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TBC1D1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human TBC1D1
Target Name	TBC1D1
Other Names	TBC; TBC1
Accession No.	Swiss-Prot#: Q86TI0Gene ID: 23216
Uniprot	Q86TI0
GeneID	23216;
Concentration	1.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43157 at dilution 1/35,

Background

TBC1D1 is the founding member of a family of proteins sharing a 180- to 200-amino acid TBC domain presumed to have a role in regulating cell growth and differentiation. These proteins share significant homology with TRE2.

Note: This product is for in vitro research use only