

Recombinant Human SEC11A

Catalog No: #GP13101



Package Size: #GP13101-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SEC11A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-179 amino acids of human SEC11A
Target Name	SEC11 homolog A, signal peptidase complex subunit
Other Names	SPC18; SPCS4A; SEC11L1; sid2895; 1810012E07Rik
Accession No.	Swissprot:P67812Gene Accession:BC014508
Uniprot	P67812
GeneID	23478;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the peptidase S26B family. The encoded protein is an 18kDa subunit of the signal peptidase complex and has been linked to cell migration and invasion, gastric cancer and lymph node metastasis. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only