

Recombinant Human NIPSNAP3A

Catalog No: #GP13121



Package Size: #GP13121-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NIPSNAP3A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-247 amino acids of human NIPSNAP3A
Target Name	nipsnap homolog 3A
Other Names	TASSC; HSPC299; NIPSNAP4
Accession No.	Swissprot:Q9UFN0Gene Accession:BC005935
Uniprot	Q9UFN0
GeneID	25934;
Storage	-20~-80°C, pH 7.6 PBS

Background

NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only