

# Recombinant Human MELTF

Catalog No: #GP13137



Package Size: #GP13137-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MELTF
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 20-302 amino acids of human MELTF
Target Name	melanotransferrin
Other Names	MTf; MFI2; MTF1; CD228; MAP97
Accession No.	Swissprot:P08582Gene Accession:BC001875
Uniprot	P08582
GeneID	4241;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a cell-surface glycoprotein found on melanoma cells. The protein shares sequence similarity and iron-binding properties with members of the transferrin superfamily. The importance of the iron binding function has not yet been identified. This gene resides in the same region of chromosome 3 as members of the transferrin superfamily. Alternative splicing results in two transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only