

# Recombinant Human SEPT1

Catalog No: #GP13141



Package Size: #GP13141-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SEPT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 168-367 amino acids of human SEPT1
Target Name	septin 1
Other Names	LARP; SEP1; DIFF6; PNUTL3
Accession No.	Swissprot:Q8WYJ6Gene Accession:BC012161
Uniprot	Q8WYJ6
GeneID	1731;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis and the maintenance of cellular morphology. This gene encodes a protein that can form homo- and heterooligomeric filaments, and may contribute to the formation of neurofibrillary tangles in Alzheimer's disease. Alternatively spliced transcript variants have been found but the full-length nature of these variants has not been determined.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only