

Recombinant Human SNAPC1

Catalog No: #GP13142



Package Size: #GP13142-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SNAPC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 169-368 amino acids of human SNAPC1
Target Name	small nuclear RNA activating complex polypeptide 1
Other Names	SNAP43; PTFgamma
Accession No.	Swissprot:Q16533Gene Accession:BC014984
Uniprot	Q16533
GeneID	6617;
Storage	-20~-80°C, pH 7.6 PBS

Background

Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only