

Recombinant Human RMND5A

Catalog No: #GP13149



Package Size: #GP13149-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RMND5A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 192-391 amino acids of human RMND5A
Target Name	required for meiotic nuclear division 5 homolog A
Other Names	CTLH; GID2; RMD5; GID2A; p44CTLH
Accession No.	Swissprot:Q9H871 Gene Accession:BC012165
Uniprot	Q9H871
GeneID	64795;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only