

Recombinant Human ADRM1

Catalog No: #GP13151



Package Size: #GP13151-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADRM1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 208-407 amino acids of human ADRM1
Target Name	adhesion regulating molecule 1
Other Names	ARM1; ARM-1; GP110
Accession No.	Swissprot:Q16186Gene Accession:BC017245
Uniprot	Q16186
GeneID	11047;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the adhesion regulating molecule 1 protein family. The encoded protein is a component of the proteasome where it acts as a ubiquitin receptor and recruits the deubiquitinating enzyme, ubiquitin carboxyl-terminal hydrolase L5. Increased levels of the encoded protein are associated with increased cell adhesion, which is likely an indirect effect of this intracellular protein. Dysregulation of this gene has been implicated in carcinogenesis. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only