

Recombinant Human TTLL1

Catalog No: #GP13155



Package Size: #GP13155-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TTLL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 224-423 amino acids of human TTLL1
Target Name	tubulin tyrosine ligase like 1
Other Names	C22orf7; HS323M22B
Accession No.	Swissprot:O95922Gene Accession:BC014968
Uniprot	O95922
GeneID	25809;
Storage	-20~-80°C, pH 7.6 PBS

Background

Catalytic subunit of the neuronal tubulin polyglutamylase complex. Modifies alpha- and beta-tubulin, generating side chains of glutamate on the gamma-carboxyl groups of specific glutamate residues within the C-terminal tail of alpha- and beta-tubulin (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only