

## Recombinant Human SMOC1

Catalog No: #GP13159



Package Size: #GP13159-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SMOC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 236-435 amino acids of human SMOC1
Target Name	SPARC related modular calcium binding 1
Other Names	OAS
Accession No.	Swissprot:Q9H4F8Gene Accession:BC011548
Uniprot	Q9H4F8
GeneID	64093;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a multi-domain secreted protein that may have a critical role in ocular and limb development. Mutations in this gene are associated with microphthalmia and limb anomalies. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only