

Recombinant Human PWP1

Catalog No: #GP13167



Package Size: #GP13167-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PWP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 302-501 amino acids of human PWP1
Target Name	PWP1 homolog, endonuclein
Other Names	IEF-SSP-9502
Accession No.	Swissprot:Q13610Gene Accession:BC001652
Uniprot	Q13610
GeneID	11137;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene contains several WD-40 repeats and is found mostly in the nucleus. The expression and localization of this protein are cell cycle dependent. Expression of this gene is upregulated in pancreatic adenocarcinoma. Three transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only