

## Recombinant Human RABGGTA

Catalog No: #GP13176



Package Size: #GP13176-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RABGGTA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 368-567 amino acids of human RABGGTA
Target Name	Rab geranylgeranyltransferase alpha subunit
Other Names	PTAR3
Accession No.	Swissprot:Q92696Gene Accession:BC003093
Uniprot	Q92696
GeneID	5875;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Catalyzes the transfer of a geranylgeranyl moiety from geranylgeranyl diphosphate to both cysteines of Rab proteins with the C-terminal sequence -XXCC, -XCXC and -CCXX, such as RAB1A, RAB3A, RAB5A and RAB7A.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only