

Recombinant Human NUP58

Catalog No: #GP13178



Package Size: #GP13178-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NUP58
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 400-599 amino acids of human NUP58
Target Name	nucleoporin 58
Other Names	NUP45; NUPL1; PRO2463
Accession No.	Swissprot:Q9BVL2Gene Accession:BC001104
Uniprot	Q9BVL2
GeneID	9818;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the nucleoporin family that shares 87% sequence identity with rat nucleoporin p58. The protein is localized to the nuclear rim and is a component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only