

# Recombinant Human ATMIN

Catalog No: #GP13184



Package Size: #GP13184-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ATMIN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 624-823 amino acids of human ATMIN
Target Name	ATM interactor
Other Names	ASCIZ; ZNF822
Accession No.	Swissprot:O43313Gene Accession:BC002701
Uniprot	O43313
GeneID	23300;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Transcription factor. Plays a crucial role in cell survival and RAD51 foci formation in response to methylating DNA damage. Involved in regulating the activity of ATM in the absence of DNA damage. May play a role in stabilizing ATM. Binds to the DYNLL1 promoter and activates its transcription.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only