

Recombinant Human PRPF3

Catalog No: #GP13185



Package Size: #GP13185-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PRPF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 484-683 amino acids of human PRPF3
Target Name	pre-mRNA processing factor 3
Other Names	PRP3; RP18; HPRP3; Prp3p; HPRP3P; SNRNP90
Accession No.	Swissprot:O43395Gene Accession:BC000184
Uniprot	O43395
GeneID	9129;
Storage	-20~-80°C, pH 7.6 PBS

Background

The removal of introns from nuclear pre-mRNAs occurs on complexes called spliceosomes, which are made up of 4 small nuclear ribonucleoprotein (snRNP) particles and an undefined number of transiently associated splicing factors. This gene product is one of several proteins that associate with U4 and U6 snRNPs. Mutations in this gene are associated with retinitis pigmentosa-18.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only