

Recombinant Human RBM15

Catalog No: #GP13200



Package Size: #GP13200-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RBM15
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 778-977 amino acids of human RBM15
Target Name	RNA binding motif protein 15
Other Names	OTT; OTT1; SPEN
Accession No.	Swissprot:Q96T37Gene Accession:BC006397
Uniprot	Q96T37
GeneID	64783;
Storage	-20~-80°C, pH 7.6 PBS

Background

Members of the SPEN (Split-end) family of proteins, including RBM15, have repressor function in several signaling pathways and may bind to RNA through interaction with spliceosome components (Hiriart et al., 2005 [PubMed 16129689]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only