

# Recombinant Human RBM19

Catalog No: #GP13203



Package Size: #GP13203-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RBM19
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 761-960 amino acids of human RBM19
Target Name	RNA binding motif protein 19
Other Names	Mrd1
Accession No.	Swissprot:Q9Y4C8Gene Accession:BC004289
Uniprot	Q9Y4C8
GeneID	9904;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only