

# Recombinant Human NID2

Catalog No: #GP13204



Package Size: #GP13204-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NID2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1237-1375 amino acids of human NID2
Target Name	nidogen 2
Other Names	NID-2
Accession No.	Swissprot:Q14112Gene Accession:BC035608
Uniprot	Q14112
GeneID	22795;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only