

Recombinant Human IPO11

Catalog No: #GP13206



Package Size: #GP13206-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IPO11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 776-975 amino acids of human IPO11
Target Name	importin 11
Other Names	RanBP11
Accession No.	Swissprot:Q9UI26Gene Accession:BC033776
Uniprot	Q9UI26
GeneID	51194;
Storage	-20~-80°C, pH 7.6 PBS

Background

Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only