

Recombinant Human EMC1

Catalog No: #GP13209



Package Size: #GP13209-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EMC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-221 amino acids of human EMC1
Target Name	ER membrane protein complex subunit 1
Other Names	CAVIPMR; KIAA0090
Accession No.	Swissprot:Q8N766Gene Accession:BC034589
Uniprot	Q8N766
GeneID	23065;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a single-pass type I transmembrane protein, which is a subunit of the endoplasmic reticulum membrane protein complex (EMC). Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only