

# Recombinant Human PLCL2

Catalog No: #GP13210



Package Size: #GP13210-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PLCL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 928-1127 amino acids of human PLCL2
Target Name	phospholipase C like 2
Other Names	PLCE2
Accession No.	Swissprot:Q9UPR0Gene Accession:BC036392
Uniprot	Q9UPR0
GeneID	23228;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play an role in the regulation of Ins(1,4,5)P3 around the endoplasmic reticulum.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only