

# Recombinant Human ZNF471

Catalog No: #GP13213



Package Size: #GP13213-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ZNF471
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 427-626 amino acids of human ZNF471
Target Name	zinc finger protein 471
Other Names	ERP1; Z1971
Accession No.	Swissprot:Q9BX82Gene Accession:BC125221
Uniprot	Q9BX82
GeneID	57573;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May be involved in transcriptional regulation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only