

# Recombinant Human SNRPD1

Catalog No: #GP13218



Package Size: #GP13218-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SNRPD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-119 amino acids of human SNRPD1
Target Name	small nuclear ribonucleoprotein D1 polypeptide
Other Names	SMD1; SNRPD; Sm-D1; HsT2456
Accession No.	Swissprot:P62314Gene Accession:BC001721
Uniprot	P62314
GeneID	6632;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a small nuclear ribonucleoprotein that belongs to the SNRNP core protein family. The protein may act as a charged protein scaffold to promote SNRNP assembly or strengthen SNRNP-SNRNP interactions through nonspecific electrostatic contacts with RNA. Two transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only