

# Recombinant Human NUTF2

Catalog No: #GP13223



Package Size: #GP13223-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NUTF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-127 amino acids of human NUTF2
Target Name	nuclear transport factor 2
Other Names	NTF2; PP15; NTF-2
Accession No.	Swissprot:P61970Gene Accession:BC002348
Uniprot	P61970
GeneID	10204;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only