

# Recombinant Human CPLX1

Catalog No: #GP13226



Package Size: #GP13226-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CPLX1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-134 amino acids of human CPLX1
Target Name	complexin 1
Other Names	CPX1; CPX-I
Accession No.	Swissprot:O14810Gene Accession:BC002471
Uniprot	O14810
GeneID	10815;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Proteins encoded by the complexin/synaphin gene family are cytosolic proteins that function in synaptic vesicle exocytosis. These proteins bind syntaxin, part of the SNAP receptor. The protein product of this gene binds to the SNAP receptor complex and disrupts it, allowing transmitter release.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only