

## Recombinant Human GOLGA7

Catalog No: #GP13229



Package Size: #GP13229-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GOLGA7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-137 amino acids of human GOLGA7
Target Name	golgin A7
Other Names	GCP16; GOLGA7A; HSPC041; GOLGA3AP1
Accession No.	Swissprot:Q7Z5G4Gene Accession:BC012032
Uniprot	Q7Z5G4
GeneID	51125;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May be involved in protein transport from Golgi to cell surface. The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only