

Recombinant Human LMO3

Catalog No: #GP13232



Package Size: #GP13232-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LMO3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-145 amino acids of human LMO3
Target Name	LIM domain only 3
Other Names	RBTN3; RHOM3; RBTN2; Rhom-3
Accession No.	Swissprot:Q8TAP4Gene Accession:BC026311
Uniprot	Q8TAP4
GeneID	55885;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the rhombotin family of cysteine-rich LIM domain oncogenes. This gene is predominantly expressed in the brain. Related family members, LMO1 and LMO2 on chromosome 11, have been reported to be involved in chromosomal translocations in T-cell leukemia. Many alternatively spliced transcript variants have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only