

Recombinant Human MTHFD2L

Catalog No: #GP13234



Package Size: #GP13234-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MTHFD2L
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 136-268 amino acids of human MTHFD2L
Target Name	methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 2 like
Accession No.	Swissprot:Q9H903Gene Accession:BC032771
Uniprot	Q9H903
GeneID	441024;
Storage	-20~-80°C, pH 7.6 PBS

Background

The specific function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only