

# Recombinant Human RNASE4

Catalog No: #GP13237



Package Size: #GP13237-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RNASE4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 29-147 amino acids of human RNASE4
Target Name	ribonuclease A family member 4
Other Names	RAB1; RNS4
Accession No.	Swissprot:P34096Gene Accession:BC015520
Uniprot	P34096
GeneID	6038;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene belongs to the pancreatic ribonuclease family. It plays an important role in mRNA cleavage and has marked specificity towards the 3' side of uridine nucleotides. Alternative splicing results in four transcript variants encoding the same protein. This gene and the gene that encodes angiogenin share promoters and 5' exons. Each gene splices to a unique downstream exon that contains its complete coding region.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only