

Recombinant Human NPC2

Catalog No: #GP13239



Package Size: #GP13239-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NPC2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 20-151 amino acids of human NPC2
Target Name	NPC intracellular cholesterol transporter 2
Other Names	HE1; EDDM1
Accession No.	Swissprot:P61916Gene Accession:BC002532
Uniprot	P61916
GeneID	10577;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein containing a lipid recognition domain. The encoded protein may function in regulating the transport of cholesterol through the late endosomal/lysosomal system. Mutations in this gene have been associated with Niemann-Pick disease, type C2 and frontal lobe atrophy.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only