

Recombinant Human ARPC5L

Catalog No: #GP13241



Package Size: #GP13241-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARPC5L
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-153 amino acids of human ARPC5L
Target Name	actin related protein 2/3 complex subunit 5 like
Other Names	ARC16-2
Accession No.	Swissprot:Q9BPX5Gene Accession:BC000018
Uniprot	Q9BPX5
GeneID	81873;
Storage	-20~-80°C, pH 7.6 PBS

Background

May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only