

Recombinant Human GSG1L

Catalog No: #GP13258



Package Size: #GP13258-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GSG1L
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 156-331 amino acids of human GSG1L
Target Name	GSG1 like
Other Names	PRO19651
Accession No.	Swissprot:Q6UXU4Gene Accession:BC016460
Uniprot	Q6UXU4
GeneID	146395;
Storage	-20~-80°C, pH 7.6 PBS

Background

As a component of the inner core of AMPAR complex, modifies AMPA receptor (AMPA) gating.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only