

Recombinant Human MPV17

Catalog No: #GP13259



Package Size: #GP13259-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MPV17
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 39-176 amino acids of human MPV17
Target Name	MPV17, mitochondrial inner membrane protein
Other Names	SYM1; MTDPS6
Accession No.	Swissprot:P39210Gene Accession:BC001115
Uniprot	P39210
GeneID	4358;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a mitochondrial inner membrane protein that is implicated in the metabolism of reactive oxygen species. Mutations in this gene have been associated with the hepatocerebral form of mitochondrial DNA depletion syndrome (MDDS).?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only