

Recombinant Human SMUG1

Catalog No: #GP13260



Package Size: #GP13260-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SMUG1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-135 amino acids of human SMUG1
Target Name	single-strand-selective monofunctional uracil-DNA glycosylase 1
Other Names	FDG; UNG3; HMUDG
Accession No.	Swissprot:Q53HV7Gene Accession:BC000417
Uniprot	Q53HV7
GeneID	23583;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that participates in base excision repair by removing uracil from single- and double-stranded DNA. Many alternatively spliced transcript variants exist for this gene; the full-length nature is known for some but not all of the variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only