

Recombinant Human FMNL2

Catalog No: #GP13261



Package Size: #GP13261-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FMNL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 880-1059 amino acids of human FMNL2
Target Name	formin like 2
Other Names	FHOD2
Accession No.	Swissprot:Q96PY5Gene Accession:BC036492
Uniprot	Q96PY5
GeneID	114793;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. Alternatively spliced transcript variants encoding different isoforms have been described but their full-length nature has yet to be determined.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only