

Recombinant Human PDAP1

Catalog No: #GP13264



Package Size: #GP13264-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDAP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-181 amino acids of human PDAP1
Target Name	PDGFA associated protein 1
Other Names	PAP; PAP1; HASPP28
Accession No.	Swissprot:Q13442Gene Accession:BC000684
Uniprot	Q13442
GeneID	11333;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a phosphoprotein that may upregulate the PDGFA-stimulated growth of fibroblasts and also downregulate the mitogenicity of PDGFB. The encoded protein in rodents has been shown to bind PDGFA with a low affinity.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only