

Recombinant Human TMEM9

Catalog No: #GP13266



Package Size: #GP13266-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TMEM9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 21-183 amino acids of human TMEM9
Target Name	transmembrane protein 9
Other Names	DERM4; TMEM9A
Accession No.	Swissprot:Q9P0T7Gene Accession:BC001106
Uniprot	Q9P0T7
GeneID	252839;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in intracellular transport.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only