

Recombinant Human IFT27

Catalog No: #GP13270

Package Size: #GP13270-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IFT27
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-186 amino acids of human IFT27
Target Name	intraflagellar transport 27
Other Names	RAYL; BBS19; RABL4
Accession No.	Swissprot:Q9BW83Gene Accession:BC000566
Uniprot	Q9BW83
GeneID	11020;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a GTP-binding protein that is a core component of the intraflagellar transport complex B. Characterization of the similar Chlamydomonas protein indicates a function in cell cycle control. Alternative splicing of this gene results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only