

Recombinant Human EEF1AKMT2

Catalog No: #GP13273



Package Size: #GP13273-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EEF1AKMT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-187 amino acids of human EEF1AKMT2
Target Name	EEF1A lysine methyltransferase 2
Other Names	Efm4; METTL10; C10orf138
Accession No.	Swissprot:Q5JPI9Gene Accession:BC026167
Uniprot	Q5JPI9
GeneID	399818;
Storage	-20~-80°C, pH 7.6 PBS

Background

Protein-lysine methyltransferase that selectively catalyzes the trimethylation of EEF1A at 'Lys-318'.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only