

Recombinant Human ZNF581

Catalog No: #GP13276

Package Size: #GP13276-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZNF581
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-197 amino acids of human ZNF581
Target Name	zinc finger protein 581
Other Names	HSPC189
Accession No.	Swissprot:Q9P0T4Gene Accession:BC010097
Uniprot	Q9P0T4
GeneID	51545;
Storage	-20~-80°C, pH 7.6 PBS

Background

ZNF581 contains 4 C2H2-type zinc fingers. ZNF581 may be involved in transcriptional regulation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only