

# Recombinant Human ASB6

Catalog No: #GP13277

Package Size: #GP13277-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ASB6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-134 amino acids of human ASB6
Target Name	ankyrin repeat and SOCS box containing 6
Accession No.	Swissprot:Q9NWX5Gene Accession:BC001719
Uniprot	Q9NWX5
GeneID	140459;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene belongs to a family of ankyrin repeat proteins that, along with four other protein families, contain a C-terminal SOCS box motif. Growing evidence suggests that the SOCS box, similar to the F-box, acts as a bridge between specific substrate-binding domains and the more generic proteins that comprise a large family of E3 ubiquitin protein ligases. Alternatively spliced transcript variants have been identified for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only